



Restriction Map and Multiple Cloning Site (MCS) of pIRES2-DsRed-Express Vector. Unique restriction sites are in bold.

Description

pIRES2-DsRed-Express contains the internal ribosome entry site (IRES; 1, 2) of the encephalomyocarditis virus (ECMV) between the MCS and a mutant of the *Discosoma* sp. red fluorescent protein DsRed-Express coding region (3–5). This permits both the gene of interest (cloned into the MCS) and the DsRed-Express gene to be translated from a single bicistronic mRNA. pIRES2-DsRed-Express is designed for the efficient selection by flow cytometry of transiently transfected mammalian cells expressing DsRed-Express and the protein of interest. This vector can also be used to express DsRed-Express alone or to obtain stably transfected cell lines without time-consuming drug and clonal selection.

DsRed-Express is a human codon-optimized (4) DsRed variant of *Discosoma* sp. red fluorescent protein (6). Its coding sequence contains nine amino acid substitutions which improve solubility of the protein and reduce the time from transfection to detection of red fluorescence (excitation and emission maxima = 557 nm and 579 nm, respectively). In addition, these substitutions reduce the level of residual green emission (5). Although DsRed-Express most likely forms the same tetrameric structure as wild-type DsRed, DsRed-Express displays a reduced tendency to aggregate (3, 5).

The MCS in pIRES2-DsRed-Express is between the immediate early promoter of cytomegalovirus ($P_{CMV IE}$) and the IRES sequence. SV40 polyadenylation signals downstream of the DsRed-Express gene direct proper processing of the 3' end of the bicistronic mRNA. The vector backbone also contains an SV40 origin for replication in mammalian cells expressing the SV40 T antigen. A neomycin-resistance cassette (Neo^r), consisting of the SV40 early promoter, the neomycin/kanamycin resistance gene of Tn5, and polyadenylation signals from the herpes simplex virus thymidine kinase (HSV TK) gene, allows stably transfected eukaryotic cells to be selected using G418. A bacterial promoter upstream of this cassette expresses kanamycin resistance in *E. coli*. The pIRES2-DsRed-Express backbone also provides a pUC origin of replication for propagation in *E. coli* and an f1 origin for single-stranded DNA production.

Use

pIRES2-DsRed-Express can be used to quickly identify cells expressing a gene of interest by screening for DsRed-Express fluorescence. Genes inserted into the MCS should include the initiating ATG codon. When DsRed-Express is expressed in mammalian cell cultures, red-emitting cells can be detected by flow cytometry or microscopy 24 hours after transfection. However, in some cases, up to 48 hours may be required for detection of red-emitting cells. pIRES2-DsRed-Express and its derivatives can be introduced into mammalian cells using any standard transfection method. If required, stable transformants can be selected using G418 (7).

Please refer to the BD Living Colors™ User Manual Volume II (PT3404-1) provided with this vector for additional information on detection of DsRed-Express.

Location of features

- Human cytomegalovirus (CMV) immediate early promoter: 1–589
Enhancer region: 59–465; TATA box: 554–560; Transcription start point: 583
C→G mutation to remove *Sac* I site: 569
- MCS: 591–665
- IRES sequence: 666–1250
- *Discosoma sp.* Red Fluorescent Protein (DsRed-Express) gene
Start codon (ATG): 1254–1256; Stop codon: 1929–1931
- SV40 early mRNA polyadenylation signal
Polyadenylation signals: 2083–2088 & 2112–2117; mRNA 3' ends: 2121 & 2133
- f1 single-strand DNA origin: 2180–2635 (Packages the noncoding strand of DsRed-Express.)
- Bacterial promoter for expression of *Kan^r* gene:
–35 region: 2697–2702; –10 region: 2720–2725
Transcription start point: 2732
- SV40 origin of replication: 2976–3111
- SV40 early promoter/enhancer
72-bp tandem repeats: 2809–2952; 21-bp repeats (3): 2956–3019
Early promoter element: 3032–3038
- Kanamycin/neomycin resistance gene: 3160–3954
G→A mutation to remove *Pst* I site: 3342; C→A (Arg to Ser) mutation to remove *Bss*H II site: 3688
- Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signals: 4190–4208
- pUC plasmid replication origin: 4539–5182

Propagation in *E. coli*

- Suitable host strains: DH5 α , HB101, and other general purpose strains. Single-stranded DNA production requires a host containing an F plasmid such as JM101 or XL1-Blue.
- Selectable marker: plasmid confers resistance to kanamycin (50 μ g/ml) to *E. coli* hosts.
- *E. coli* replication origin: pUC
- Copy number: ~500
- Plasmid incompatibility group: pMB1/ColE1

References

1. Jackson, R. J., et al. (1990) *Trends Biochem. Sci.* **15**:477–483.
2. Jang, S. K., et al. (1990) *J. Virol.* **62**:2636–2643.
3. BD Living Colors DsRed-Express (July 2002) *Clontechniques XVII*(3):16–17.
4. Matz, M. V., et al. (1999) *Nature Biotechnol.* **17**:969–973.
5. Bevis, B. J. & Glick, B. S. (2002) *Nature Biotechnol.* **20**:83–87.
6. Haas, J., et al. (1996) *Curr. Biol.* **6**:315–324.
7. Gorman, C. (1985). In *DNA cloning: A practical approach*, vol. II. Ed. D.M. Glover. (IRL Press, Oxford, U.K.) pp. 143–190.

Notice to Purchaser

This product is the subject of pending U.S. patents.

Not-For-Profit-Entities: Orders may be placed in the normal manner by contacting your local representative or BD Biosciences Clontech Customer Service at either 800-662-2566 or 650-424-8222, extension 1. At its discretion, BD Biosciences Clontech grants not-for-profit research entities a worldwide, non-exclusive, royalty-free, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product or its components to third parties. Any other use of this product will require a license from BD Biosciences Clontech. Please contact the Product Manager for Cell Biology at either 800-662-2566 or 650-424-8222, extension 7816.

For-Profit entities that wish to use this product in non-commercial or commercial applications are required to obtain a license from BD Biosciences Clontech. For license information please contact the Product Manager for Cell Biology at either 800-662-2566 or 650-424-8222, extension 7816.

Use of the IRES sequence is covered by U.S. Patent No. 4,937,190 and is limited to use solely for research purposes. Any other use of the IRES sequence requires a license from Wisconsin Alumni Research Foundation.

The CMV promoter is covered under U.S. Patent Nos. 5,168,062, and 5,385,839 assigned to the University of Iowa Research Foundation.

The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by BD Biosciences Clontech. This vector has not been completely sequenced.

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. BD Biosciences Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of BD Biosciences Clontech.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2004 BD