



```

323
•
GAATTCGAGCTCGGTACCCGGGGATCCTCTAGTCAGCTGACGCGT
EcoRI SacI KpnI BamHI PvuII MluI
368
•
GCTAGCGCGGCCGCATCGATAAGCTTGTGCGACGATATCTCTAGA
NheI EagI ClaI HindIII SalI EcoRV XbaI
NotI AccI

```

Map and Multiple Cloning Site (MCS) of pTRE-Tight Vector. Unique restriction sites are in bold.

Description

pTRE-Tight is a response plasmid that can be used to express a gene of interest (Gene X) in our BD™ Tet-On and BD™ Tet-Off Gene Expression Systems and Cell Lines (1). The Tet Expression Systems and Cell Lines give researchers ready access to the tetracycline-regulated expression systems described by Gossen & Bujard (2; Tet-Off) and Gossen *et al.* (3; Tet-On). pTRE-Tight contains an MCS immediately downstream of the **Tet-responsive P_{tight} promoter**. cDNAs or genes inserted into the MCS will be responsive to the **tTA and rtTA regulatory proteins** in the Tet-Off and Tet-On systems, respectively. P_{tight} contains a modified Tet response element (TRE_{mod}), which consists of seven direct repeats of a 36-bp sequence that contains the 19-bp tet operator sequence (*tetO*). The TRE_{mod} is just upstream of the minimal CMV promoter ($P_{minCMV\Delta}$), **which lacks the enhancer that is part of the complete CMV promoter**. Consequently, P_{tight} is silent in the absence of binding of TetR or rTetR to the *tetO* sequences. Note that the cloned insert must have an initiating ATG codon. In some cases, addition of a **Kozak consensus ribosome binding site** (4) may improve expression levels; however, many cDNAs have been efficiently expressed in Tet systems without the addition of a Kozak sequence. pTRE-Tight-Gene X plasmids should be cotransfected with the Linear Hygromycin Marker (#631625 or #6202-1, not included) or Linear Puromycin Marker (#631626 or #6203-1, not included) to permit selection of stable transfectants. pTRE-Tight was derived from pTRE, originally described as pUHD10-3 in reference 5.

The pTRE-Tight-Luc Control Vector, packaged with the pTRE-Tight Vector, contains an additional 1,649 bp encoding firefly luciferase inserted into the MCS. This vector can be used as a reporter of induction efficiency using standard luciferase detection reagents. It is not intended as a cloning vector.

Location of features

- P_{tight} Tet-responsive promoter: 3–318
 - Tet response element (TRE_{mod}): 3–252
 - Location of seven *tetO* 19-mers: 12–30; 48–66; 83–101; 119–137; 155–173; 190–208 & 226–244
 - Fragment containing $P_{\text{minCMV}\Delta}$: 258–317
 - TATAA box: 280–286
- Multiple cloning site (MCS): 323–411
- Fragment containing SV40 poly A signal: 406–606
- Fragment containing Col E1 origin of replication: 780–1379
- Ampicillin resistance gene (β -lactamase): 2536–1541

Propagation in *E. coli*

- Suitable host strains: DH5 α and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (100 $\mu\text{g/ml}$) in *E. coli* hosts.
- *E. coli* replication origin: Col E1

References

1. HT1080 Cell Line & pTRE2 Vector (January 1999) *Clontechniques* **XIV**(1):23.
2. Gossen, M. & Bujard, H. (1992) *Proc. Natl. Acad. Sci USA* **89**:5547–5551.
3. Gossen, M., *et al.* (1995) *Science* **268**:1766–1769.
4. Kozak, M. (1987) *Nucleic Acids Res.* **15**:8125–8148.
5. Resnitzky, D., *et al.* (1994) *Mol. Cell. Biol.* **14**:1669–1679.

Note:

The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by BD Biosciences Clontech. This vector has not been completely sequenced.

Notice to Purchaser

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. BD Biosciences Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of BD Biosciences Clontech.

Use of the Tetracycline controllable expression systems (the "Tet Technology") is covered by a series of patents including U.S. patents #5,464,758 and #5,814,618 which are proprietary to Abbott Laboratories. Academic research institutions are granted an automatic license with the purchase of this product to use the Tet Technology only for internal, academic research purposes, which license specifically excludes the right to sell, or otherwise transfer, the Tet Technology or its component parts to third parties. In accepting this license, all users acknowledge that the Tet Technology is experimental in nature. Abbott makes no warranties, express or implied or of any kind, and hereby disclaims any warranties, representations, or guarantees of any kinds as to the Tet Technology, patents, or products. All others are invited to request a license from Abbott prior to purchasing these reagents or using them for any purpose. BD Biosciences Clontech is required by its licensing agreement to submit a report of all purchasers of the Tet-controllable expression systems to Abbott. For license information, please contact:

Abbott Bioresearch Center
100 Research Drive,
Worcester, MA 01605-4314, U.S.A.,
Fax: 508-755-8506

BD, BD logo and all other trademarks are the property of Becton, Dickinson and Company.

© 2003 BD